

# PROVOX<sup>®</sup>

## GuideWire

**MD**

**Rx  
ONLY**

**CE** 0413

**STERILE | EO**

42°C  
108°F  
RT  
2°C  
36°F

Fig.1


Fig. 2







Fig. 3










Fig. 4

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and trademarks**

Provox® is a registered trademark owned by Atos Medical AB, Sweden.

For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Produsent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici; მწარმოებელი; Ἀρτοποιός; İstehsalçı; Производитель; Produzen; Pengilang; 製造業者; 제조사; 製造商; 制造商; הגורם; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productiedatum; Date de fabrication; Data di fabbricazione; Fecha de fabricación; Data de fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato; Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev; Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία κατασκευής; Дата на производство; Üretim tarihi; წარმოების თარიღი; Ἀρτοποιημάτων ἡμερομηνία; İstehsal tarihi; Дата изготовления; Tanggal produksi; Tarikh pembuatan; 製造日; 제조일자; 製造日期; 製造日期; תאריך ייצור; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite d'utilisation; Data di scadenza; Fecha de caducidad; Data de validade; Används senast; Anvendes inden; Brukes innen; Viimeinen käyttöpäivä; Sídasti notkunardagur; Kõlblik kuni; Derīguma termiņš; Naudoti iki; Datum expirace; Felhasználatóság dátuma; Dátum najneskoršej spotreby; Rok uporabnosti; Data ważności; Data expirării; Rok uporabe; Rok upotrebe; Ημερομηνία λήξης; Срок на годност; Son kullanım tarihi; ვარგისობა; ἡμερομηνία ἕως ἧς; Son istifadə tarixi; Срок годности; Gunakan sebelum tanggal; Tarikh guna sebelum; 使用期限; 사용 기한; 使用截止日期; 使用截止日期; תאריך אחרון לשימוש; تاريخ انتهاء الصلاحية

Batch code; Chargencode; Partijnummer; Code de lot; Codice lotto; Código de lote; Código do lote; Batchkod; Batchkode; Batchkode; Eränumero; Lotunúmer; Partii kood; Partijas kods; Partijos kodas; Kód šarže; Sarzs kód; Kód šarže; Številka serije; Kod partii; Codul lotului; Šifra serije; Šifra serije; Κωδικός παρτίδας; Код на партида; Parti kodu; სერის კოდი; Isürpawarınalıq kodu; Partiya kodu; Код партии; Kode batch; Kod kelompok; バッチコード; 배치 코드; 批次代碼; 批号; קוד אצווה; كود الدفعة

Product reference number; Produktreferenznummer; Referentienummer product; Numéro de référence du produit; Numero riferimento prodotto; Número de referencia del producto; Referência do produto; Referensnummer; Produktreferenznummer; Produktreferenznummer; Artikkelnummer; Tilvísunarnúmer vöru; Toote viitenumber; Produkta atsauces numurs; Gaminio katalogo numeris; Referenční číslo produktu; Termék referenciaszáma; Referenčné číslo produktu; Referenčna številka izdelka; Numer katalogowy produktu; Numărul de referință al produsului; Referentni broj proizvoda; Referentni broj proizvoda; Αριθμός αναφοράς προϊόντος; Референтен номер на продукт; Ürün referans numarası; პროდუქტის საკონტროლო ნომერი; Ἀρτοποιημάτων ἡμερομηνία; Məhsula istinad nömrəsi; Справочный номер изделия; Nomor referensi produk; Nombor rujukan produk; 製品参照番号; 제품 참조 번호; 產品參考號; 产品参考编号; מספר סימוכין של המוצר; الرقم المرجعي للمنتج

Do not re-use; Nur zum Einmalgebrauch; Niet hergebruiken; Ne pas réutiliser; Non riutilizzare; No reutilizar; Não reutilizar; Återanvänd ej; Må ikke genanvendes; Skal ikke gjenbrukes; Ei uudelleenkäytettävä; Má ekki endurnýta; Mitte taaskasutada; Neizmantot atkārtoti; Pakartotinai nenaudoti; Nepoužívejte opakovaně; Ne használja fel újra; Nepoužívajte opakovane; Ne uporabljajte ponovno; Nie używać ponownie; A nu se reutiliza; Nemojte ponovno upotrebljavati; Ne koristiti ponovo; Μην επαναχρησιμοποιείτε; Не използвайте повторно; Tekrar kullanmayın; არ გამოიყენოთ ხელახლა; Մեկանգամյա օգտագործման համար; Yeniden istifadə etməyin; Не подлежит повторному использованию; Jangan

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# 1. Descriptive Information

## 1.1 Indications for use

The **Provox GuideWire** is a sterile single use insertion device intended for placement of a sterile Provox indwelling Voice Prosthesis after total laryngectomy (primary or secondary puncture), or for retrograde replacement of a sterile Provox indwelling Voice Prosthesis.

## 1.2 Contraindications

Do not use if the patient has anatomical abnormalities, e. g. significant pharyngeal stenosis above the puncture site or severe trismus. Significant pharyngeal stenosis may preclude insertion of the voice prosthesis. Severe trismus may preclude proper protection of the pharyngeal wall during secondary puncture leading to harm of the esophageal tissue.

## 1.3 Description of the device

The Provox GuideWire is a sterile, single use device for introduction and replacement of sterile Provox indwelling Voice Prostheses. The GuideWire has a connector for attachment of the safety strap of the new voice prosthesis and an 8 mm Stopper for transoral removal of the remnant of the old voice prosthesis.

The Provox GuideWire package contains the following items:

- 1 Provox GuideWire, sterile (Fig.1),
- 1 Instructions for use Provox GuideWire, non-sterile.

The Provox GuideWire has two functional parts (Fig.1)

1. Stopper – for transoral removal of the old Provox Voice Prosthesis (made of ABS plastic) (Fig.1.1),
2. Connector – for attachment of a Provox Voice Prosthesis safety strap (made of nylon plastic material (Polyamide) (Fig.1.2).

The Provox GuideWire tube is made of a PVC plastic material.

## 1.4 WARNINGS

### Pre-surgery

- **DO NOT** use the product if the package or product is damaged or opened. Unsterile product may cause infection.
- **DO NOT REUSE** and **DO NOT RESTERILIZE** by any method. Intended for single use only. Reuse may cause cross contamination. Cleaning and resterilization may cause structural damage.

### During insertion

- **DO** always proceed slowly and without using excessive force when pulling the GuideWire through. Tissue damage or dislodgement or extrusion of parts of the Provox GuideWire may otherwise occur which can lead to subsequent ingestion or tissue damage.
- **DO NOT** proceed if the voice prosthesis gets stuck in pharynx or esophagus, this may lead to tissue damage.

For further information about these events, see section Adverse events.

## 1.5 PRECAUTIONS

Always assess the suitability before using the Provox GuideWire in the patient. In cases with lacking suitability, do not use the product.

- **DO** carefully assess patients with bleeding disorders or patients undergoing anticoagulant treatment for the risk of bleeding or hemorrhage prior to placement or replacement of the voice prosthesis.
- **DO** ensure that in case of known or suspected pharyngeal stenosis, dilatation of the pharynx should be performed before using the GuideWire.
- **DO** always use aseptic technique during the procedure in order to reduce the risk of infection.
- **DO** ensure that the safety strap of the voice prosthesis is adequately attached to the GuideWire. If proper attachment is not achieved, the voice prosthesis may come loose from the GuideWire, causing failure to complete the procedure.
- **DO NOT** load or reload the voice prosthesis if the safety strap has been cut off or is broken.

## 1.6 Adverse events and troubleshooting information

**Aspiration of parts of the Provox GuideWire or parts of the Provox System** – Accidental aspiration of parts of the Provox GuideWire or parts of the Provox System may occur. As with any other foreign body, complications from aspiration of a component may cause obstruction or infection. Immediate symptoms may include coughing, wheezing or other abnormal breathing sounds, dyspnea, and respiratory arrest, partial or inadequate air exchange and/or asymmetrical chest movement with respiration. Complications may include pneumonia, atelectasis, bronchitis, lung abscess, bronchopulmonary fistula and asthma. If the patient can breathe, coughing may remove the foreign body. Partial airway obstruction or complete airway obstruction requires immediate intervention for removal of the object. If aspiration of the object is suspected, the object should be located and retrieved endoscopically using a non-toothed grasping forceps.

**Ingestion of parts of the Provox GuideWire or parts of the Provox System** – Accidental ingestion of parts of the Provox GuideWire or parts of the Provox System may occur. As with any other foreign body, the symptoms caused by ingestion of parts of the GuideWire or parts of the Provox System, depend largely on size, location, degree of obstruction (if any) and the length of time it has been present. Ingested components that have remained in the esophagus may be removed by esophagoscopy or observed for a short period of time. The object may pass spontaneously into the stomach. Foreign bodies that pass into the stomach usually pass through the intestinal tract. Surgical removal of foreign bodies from the intestinal tract must be considered when bowel obstruction occurs, bleeding is present, perforation occurs or the object fails to pass through the intestinal tract.

Spontaneous passage of the object may be awaited for 4-6 days. The patient should be instructed to observe the stools for the ingested object. If the object does not pass spontaneously, or if there are signs of obstruction (fever, vomiting, abdominal pain) a gastroenterologist should be consulted. The device may be retrieved by using a non toothed grasping forceps.

**The Provox GuideWire gets stuck in the pharyngeal mucosal wall** – The Provox GuideWire may get stuck in or be interfered by the pharyngeal mucosal wall. With slight pressure, the Provox GuideWire will generally bend near the tip and slide upwards towards the pharynx. Stop the procedure if the general bending does not help.

**Hemorrhage/Bleeding of the puncture** – Slight bleeding from the edges of the tracheoesophageal (TE) puncture may occur during intended use of Provox GuideWire and generally resolves spontaneously. Patients on anti-coagulant therapy, however, should be carefully evaluated for the risk of hemorrhage prior to secondary puncture and prosthesis placement.

**Rupture of the tracheoesophageal tissue** – In case of a rupture of the TE tissue the TE puncture procedure should be abandoned and the rupture should be sutured immediately. TE puncture should only be repeated after proper healing of the tissues.

**The Provox GuideWire gets stuck inside the esophagus** – If the Provox GuideWire gets stuck inside the esophagus, do not use excessive force. Attempts to pull out a stuck GuideWire may cause tissue damage and/or breaking the GuideWire. The device has to be retrieved with an endoscope.

## 2. Instructions for use

### 2.1 Preparation

Check the integrity of the sterile package. Do not use the product if the package is damaged or opened.

### 2.2 Operating instruction

#### 2.2.1 Primary puncture placement (Fig 2):

1. Insert the Provox Pharynx Protector in the open pharynx (Fig 2.1).

2. Use the Provox Trocar with Cannula to create a TE-puncture (about 8-10 mm from the edge of the tracheostoma) (Fig 2.2).

3. Remove the Trocar after the puncture, while leaving the Trocar Cannula in place (Fig 2.3).

4. Introduce the Provox GuideWire in anteroposterior direction through the opening of the Trocar Cannula.

The connector of the GuideWire appears out of the upper opening of the Pharynx Protector, which then can be removed (Fig 2.4).

**CAUTION:** Do not attach the voice prosthesis when the Pharynx Protector is in place, since the voice prosthesis can get stuck.

5. The safety strap of the voice prosthesis of the proper model, diameter and size is attached and secured to the GuideWire connector. The Trocar Cannula is removed from the puncture opening and is pulled back.

Careful traction on the GuideWire will bring the strap of the prosthesis into the TE-puncture. The voice prosthesis is carefully rotated into position with two non-toothed hemostats firmly holding the tracheal flange (Fig 2.5).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

6. When the voice prosthesis is correctly positioned, the safety strap is cut off while still attached to the GuideWire and the voice prosthesis is turned so that the part where the safety strap was attached points downwards into the trachea. The laryngectomy procedure can then be completed (Fig 2.6).

#### 2.2.2 Secondary puncture placement (Fig 3):

1. Insert a rigid endoscope (or equivalent instrument) into the esophagus (Fig 3.1).

2. Use the Provox Trocar with Cannula to create a TE-puncture (about 8-10 mm from the edge of the tracheostoma) (Fig 3.2).

3. Remove the Trocar after the puncture, while leaving the Trocar Cannula in place (Fig 3.3).

4. Introduce the Provox GuideWire through the Trocar Cannula. The GuideWire is pushed upwards through the scope until it is stopped by the Stopper at the end of the Cannula. The endoscope can then be removed (Fig 3.4).

**CAUTION:** Do not attach the voice prosthesis when the endoscope is in place, since the voice prosthesis can get stuck.

5. The safety strap of the voice prosthesis of the proper model, diameter and size is attached and secured to the GuideWire connector. The Trocar Cannula is removed from the puncture opening and is pulled back.

Careful traction on the GuideWire brings the strap of the voice prosthesis into the TE-puncture. The voice prosthesis is carefully rotated into position with two non-toothed hemostats firmly holding the tracheal flange (Fig 3.5).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

6. When the voice prosthesis is correctly positioned, the safety strap is cut off while still attached to the GuideWire and the voice prosthesis is turned so that the part where the safety strap was attached points downwards into the trachea (Fig 3.6).

### 2.2.3 Replacement procedure (Fig 4):

#### **Anterograde vs. retrograde insertion**

If the sterile indwelling Voice Prosthesis cannot be inserted in an anterograde manner, retrograde insertion by use of the Provox GuideWire may be performed.

Replacement should be carried out as an outpatient procedure using local anaesthesia. Occasionally, the procedure may be carried out under general anesthesia upon determination by the physician.

#### **Anterograde removal of the old prosthesis (not Provox(1)):**

1. Grasp the tracheal flange with a non-toothed hemostat and pull the prosthesis out from the TE puncture.
2. Insert the GuideWire through the TE puncture and push the GuideWire upwards through the pharynx and out of the mouth.

#### **Retrograde removal of the old prosthesis (Only Provox(1)):**

The old Provox(1) voice prosthesis can be removed transorally by the use of the Provox GuideWire.

1. While the old voice prosthesis is still in situ, introduce the GuideWire through the old prosthesis and push the GuideWire upwards through the pharynx and out of the mouth.
2. Grasp the tracheal flange of the prosthesis with a non-toothed hemostat. Cut it off the tracheal flange with a scalpel and remove it. With the tracheal flange removed, the Stopper of the GuideWire will pull the remains of the old voice prosthesis into the esophagus and up through the mouth. Discard the old voice prosthesis.

#### **Placement of a sterile indwelling Provox Voice Prosthesis:**

The GuideWire Connector is now outside the mouth (Fig 4.1).

3. Attach and secure the new voice prosthesis to the GuideWire Connector. Pull the GuideWire downwards, towards the existing TE-puncture. Careful traction on the GuideWire brings the strap of the voice prosthesis in to the TE-puncture (Fig 4.2).

**CAUTION:** Do not pull the voice prosthesis in place with the GuideWire.

4. Introduce the tracheal flange of the voice prosthesis into the puncture with the help of a non-toothed hemostat. Then turn the voice prosthesis so that the safety strap is pointing downwards, into the trachea. Let the patient drink and swallow while observing the voice prosthesis, to ensure no leakage occurs. Cut off the safety strap while it is still attached to the GuideWire and discard the GuideWire (Fig 4.3).

## 2.3 Cleaning and sterilization

The Provox GuideWire is provided sterile (EO) and intended for single use only and can NOT be cleaned or resterilized.

## 2.4 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## 2.5 Additional devices

**Provox Trocar with Cannula:** Stainless steel instruments for performing a primary or secondary TE puncture. The Provox Trocar creates a small, rounded puncture, while the Cannula provides a guided passage for the GuideWire through the freshly made TE puncture.

**Provox Pharynx Protector:** Stainless steel instrument for protection of the pharyngeal wall during primary TE puncture.

## 3. Additional information

### 3.1 Ordering information

See back of this Instruction for Use.

### 3.2 Date of printing

See version number on the back cover of this Instruction for use.

### 3.3 User assistance information

For additional help or information please see back cover of this manual for contact information.

## 4. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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