


## Technical Data

### Tracoe Mini for neonates and infants

REF 350



| Size | ID  | OD  | A  | C  | B  | TL | OB | BA  |
|------|-----|-----|----|----|----|----|----|-----|
|      | mm  | mm  | mm | mm | mm | mm | mm | °   |
| 2.5  | 2.5 | 3.6 | 6  | 22 | 2  | 30 | 32 | 120 |
| 3.0  | 3.0 | 4.3 | 7  | 23 | 2  | 32 | 34 | 120 |
| 3.5  | 3.5 | 5.0 | 8  | 23 | 3  | 34 | 37 | 120 |
| 4.0  | 4.0 | 5.6 | 9  | 24 | 3  | 36 | 40 | 120 |

**ID:** inside diameter (clear width) of tube; **OD:** outside diameter of tube;  
**A:** length; **C:** length; **B:** length; **TL:** length along centre line from start of neck flange to bottom of tube; **OB:** length along outer bend from start of neck flange to bottom of tube; **BA:** bending angle