

Tracoe Neck Straps

Neckbands

Product Description

There are many essential products that play an important role in tracheostomy care. Our assortment of essentials includes neckbands, dressings, products for care & cleaning as well as clothing & protectors and accessories. They fulfill a wide range of user needs and preferences, are easy to use, comfortable, discreet, offer security, and are functional.

General Information

Product line	Neckbands
Legal manufacturer	Tracoe Medical GmbH Reichelsheimer Str. 1/3 55268 Nieder-Olm Germany
Country of origin	Germany
Handling	
Shelf life	5 years

Material Specifications

Material padding	PU foam
Material fabric	Cotton
Material Hook and Loop fastener	Polyamide, Nylon
Biocompatibility	ISO 10993-1, ISO 10993-5 & ISO 10993-12

Storage



Packaging

Packaging unit	10 pcs
Material outer packaging:	Resealable Plastic bag



Certification

Product classification acc. MDR	I
CE Marking	
GMDN	35752
Quality system certified to	EN ISO 13485
FDA registration	Exempt/D215988

PRODUCT DATA SHEET



Product Specifications

Product Name	REF	Length	Adjustable in length	Width	Padded	Hook and Loop fastener	
Tracoe Neck Strap One-Piece	REF 903-B-P	42 cm	No	Approx. 3.5cm	No	yes	
Tracoe Neck Strap Padded	REF 903-E-P	Max length 45 cm	yes	Approx. 3.0 cm	yes	yes	



Tracoe Medical GmbH,
Reichelsheimer Str. 1 / 3,
55268 Nieder-Olm
Germany



Breathing-Speaking-Living

Atos Medical AB
Hyllie Boulevard 17
215 32 Malmö
Sweden
www.atosmedical.com