

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	DD	Pontus Eklund - X-PONEKL	2020-04-21 - 14:58
Reviewed:	DD	Jon Berg - JONBER	2020-04-21 - 15:47
Approved:	DD	Fredrik Calais - FRECAL	2020-04-21 - 16:15
Released:	DD	Jon Berg - JONBER	2020-05-28 - 15:08

This document has been electronically signed by the persons above.

Product Information

Freevent® Neckband



Product description:

The Freevent Neckbands are tube/button holders for patients breathing through a tracheostoma. Two-piece variants connect with hook-and-loop whereas one-piece variants connect with twill ties.

Adult and Child sizes

Product Information

Document ID:	PF094-01-TechInfo	Edition:	04
Manufacturer:	Atos Medical AB Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden		
Classification: (EU) 2017/745	Class I (1.1 rule 1)		
Intended Use:	The Freevent Neckbands are used for holding a tube/button in place by wearing it around the neck and connecting the ends of the neckband to the tube/button.		
CE Mark:	Yes, the devices are CE marked.		
GMDN code:	63438 (Tracheostomy tube neck holder, single use)		
Sterilization:	Non-sterile		
Raw material:	Nylon, cotton and polyester polyurethane foam.		
Latex information:	Not manufactured with natural rubber latex		
Biological origin:	The device is not manufactured with materials derived from human or animal source.		
Handling and storage:	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.		
Waste handling and disposal:	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.		
Hazardous components:	None		
Expiration date:	5 years after manufacturing.		
Packaging:	REF 1651, 1661, 1652 and 1662: 10 pcs in a box, each neckband individually packed in a polybag (polyethylene) REF 1751, 1761, 1752 and 1762: 100 pcs in a box, each neckband individually packed in a polybag (polyethylene)		

Product Information

Devices under Basic UDI-DI: 7331791-GEN-A-000-0000-E6

REF	Name	UDI-DI
1651	Freevent Neckband, one-piece , small	07331791012914
1652	Freevent Neckband, two-piece, small	07331791012921
1661	Freevent Neckband, one-piece, large	07331791012938
1662	Freevent Neckband, two-piece, large	07331791012945
1751	Freevent Neckband, one-piece, small	07331791013607
1752	Freevent Neckband, two-piece, small	07331791013621
1761	Freevent Neckband, one-piece, large	07331791013614
1762	Freevent Neckband, two-piece, large	07331791013638



Freevent Neckband, one-piece



Freevent Neckband, two-piece

Atos Medical AB Compatible products:

Range	BASIC UDI-DI
Provox LaryTube	7331791-LTU-0-000-0002-3E
Provox LaryButton	7331791-LTU-0-000-0000-38
Provox LaryButton Sizer Kit	7331791-LTU-0-000-0001-3B